

# MedReady™



- **Timed-alarm system provides next dosage automatically**
- **Operates by easy-to-use microchip technology**
- **Holds up to 7 days of 4 times a day dosing**
- **Costs less than similar devices**
- **Flashing light available**

 **In-home Care**  
Because there's no place like home!

Automatic Medication Dispensers  
are provided free with  
Companion/Homemaker Services\*

Call Today

1-800-HomeCare

904-737-7667

or visit our website at:

[www.In-HomeCare.com](http://www.In-HomeCare.com)

Florida Lic. 30211295

\*Some Restrictions may apply

## Automated Medication Dispenser

 **In-home Care**  
Because there's no place like home!

1-800-HomeCare

904-737-7667





## FEATURES



### Basic 1600

Large 1/2" L.E.D. displays the time and aids setting alarms.

- Audible alarm can be set up to 4 times a day:
  - ✦ "Time for medicine" alarm beeps 5 seconds on, 1 second off continuously until the medication door is slid open and the "off" position is activated.
  - ✦ Note: Disabling alarm requires the medication door to be slid open.
  - ✦ Programmable alarm sounds for 30, 45, 60, or 250 minutes, as user selects
  - ✦ Powerful alarm sounds similar in tone and volume to a beeping microwave alarm but alarms on for 5 seconds and off for one second until disabled.
  - ✦ Audible alarm is 90dB, 2800-3000Hz. (1600LF is approximately 400Hz)
- Only the most recent medication dose is available. Pill tray rotates to the next dose when the next audible alarm sounds. If a patient misses a dose, the pill dispenser denies access to the missed dose.

# MedReady™

## Automated Medication Dispenser

- MedReady holds 28 doses of dry medication such as pills or capsules.
  - ✦ 4 Doses a day for 7 days
  - ✦ 3 Doses a day for 9 days
  - ✦ 2 Doses a day for 14 days
  - ✦ 1 Dose a day for 28 days
- Alarm in pill dispenser is as easy to set as a digital clock.
- Medication is under lock and key, not just a latch.
- Pill compartments are designed to fit adult fingers and "ramped" for easy removal of medication.
- Compartment divisions keep medications in place even if the unit is turned upside-down.
- MedReady should remain plugged into electric power, although there is a 48 hour rechargeable battery back-up for travel or power failures. MedReady should not be used in wet or very damp environments.
- "OFF" feature for long term disabling of the dispenser.
- "Early dosage" feature allows taking dosage with client.

## OPTIONS

for

### MedReady Medication Dispenser 1600

1. **5310 Extra Pill Tray** with clear plastic lid. Preload an additional tray for convenience.
2. **1600FL Alarm With Flashing Light** for hearing impaired. Audible alarm PLUS Flashing Light. 
3. **1600LF Lower Frequency Alarm** is louder for those losing high frequency hearing.
4. **1650 Remote Monitoring:** Simply plug the added telephone cable connector from the MedReady unit to the phone line.
  - ✦ Patient only needs phone line. Caregiver accesses secure web site to review actual medication times.
  - ✦ Monthly service plan.
  - ✦ Calls up to 3 phone numbers, leaving messages, and sends e-mail, if dose is not taken.
  - ✦ A record of the medication compliance and/or non-compliance is available on the **medready.net** web site.